

DATAGRABBER

DataGrabber, Apteryx's patented bridge to practice management programs, is designed for Windows XP, Vista, Windows 7, Server 2003, Server 2008 operating systems including Terminal Services and Citrix.

Depending on the PM application, an ODBC Data Source may need to be specified on each workstation. Contact Apteryx, Inc at **1-877-278-3799** for consultation.

Application Requirements:

- Installed Practice Management (PM) Solution
- Appropriate Apteryx DataDrill to extract patient information from PM Database
- Third Party Dental or Medical Imaging Software

i DataGrabber works with most third-party practice management software and will launch most third-party dental and medical imaging software. Contact Apteryx for the latest support information.

Minimum System Requirements

The system must conform to the minimum requirements specified by Microsoft® for the operating system installed. At minimum, DataGrabber requires Windows® XP (SP2 or above) which has the following minimum requirements:

- Pentium 233-megahertz (MHz) processor or faster (300 MHz is recommended)
- At least 64 megabytes (MB) of RAM (128 MB is recommended)
- At least 1.5 gigabytes (GB) of available space on the hard disk
- CD-ROM or DVD-ROM drive
- Keyboard and a computer mouse device, or some other compatible pointing device
- Video adapter and monitor with 1024 x 768 or higher resolution
- Sound card
- Speakers or headphones

Product Name	DataGrabber		
Department	Technical Support	Revision Date	8/1/2011
Category		Revision #	0.1
Prepared By	SN	Effective Date	-8/1/2011

This a controlled document. Printed or downloaded copies are considered uncontrolled and for reference only.